

Bathroom aids



Catálogo *online*



AYUDAS

García 1880

MOVILIDAD




www.garcia1880.com

Raised toilet seats



Toileting and daily hygiene support aids which allow you to create an accessible bathroom.

2172 Soft raised toilet seat.
Inner dimensions 21x27.



21721 Soft raised toilet seat with lock.



2173 Soft raised toilet seat with lid and lock.



2204 Adjustable toilet side safety rail. Width: 14-26 cm. Size: 62 cm length x 47 cm width x 22 cm height.



2175 Raised toilet seat with arms and lock.



2275 Raised toilet seat with flip up armrests. Height adjustable up to 15 cm.





Soft touch raised toilet seat, without lid. 150 kg.


227605 **227610** **227616**
5 cm height. 10 cm height. 15 cm height.



Soft touch raised toilet seat with lid.

227705 **227710** **227715**
5 cm height. 10 cm height. 15 cm height.

These raised toilet seats adapt easily to any kind of toilet designs.
They include locks for greater adjustment.



Raised toilet seat without lid.

227805 **227810** **227815**
5 cm height. 10 cm height. 15 cm height.



Raised toilet seat with lid.

227905 **227910** **227915**
5 cm height. 10 cm height. 15 cm height.

2280
Bidet for raised toilet seats.



22801
Portable bidet for WC.
41 x 37 x 12 cm.



395 gr.



Bathroom Grab Rails




Support products designed to improve user safety and help with mobility in a variety of tasks in the bathroom and around the home.

It is recommended that the installer values the screws to use when fixing depending on the type of wall or the place where it is installed.


2178

Steel foldable grab rail, white epoxy coated.




21781

Nickel chromed steel foldable assisting rail. 74 cm long.




21782

Stainless steel foldable assisting rail. 82 cm long.



2179

Foldable steel grab rail with height adjustable leg support.



2198015

Stainless steel assisting rail, with double U shape. 80 cm long.



2205

Portable handrail with magazine rack. Width: 78cm; height: 67, 5 cm; depth: 48, 5 cm.



2198020

Stainless Steel Grab rail wall to floor. Measures: 80x85cm.



MAX.
KG
80 kg.



70x50x12 cm



It is recommended that the installer values the screws to use when fixing depending on the type of wall or the place where it is installed.

2194031 River grab bar 31 cm.

2194041 River grab bar 41 cm.






Anti-slip grip. Height: 6,7cm.
Diameter: 32mm.




2193031 Ocean grab bar 31 cm.

2193041 Ocean grab bar 41 cm.



Corrugated and non-slip tube material:
aluminium and ABS. Diameter: 32mm.
Height: 6,7 cm.



2196031 **2196041** **2196061** **2196082**

White steel grab bars.




2197031 **2197041** **2197061** **2197082**

Chromed Steel grab bars.

2198031 31 cm stainless steel grab bar.

2198041 41 cm stainless steel grab bar.

2198061 61 cm stainless steel grab bar.





2198010
Stainless Steel angle grab bar.
Measurements: 50x70 cm.



2195061
61 cm white plastic grab bar.


2196025
White steel angle bar.

2197025
Chromed Steel angle grab bar.




Length: 30 cm.
Diameter: 32 mm.

Suction grab bar with safety indicator



2199040
40 cm corrugated plastic suction grab bar.



2199030
30 cm corrugated plastic suction grab bar.





Bath aids to remain seated during bathing or showering. They facilitate people with disabilities or reduced mobility to perform these cleaning functions.

2186

Swivel bath seat. Frame made of white epoxy coated steel.




120 kg.



5'4 kg.



73x58x32 cm



21871

Aluminium swivel seat with adjustable width.



150 kg.



56x45x26 cm

Quick assembly, no tools needed.



2183

Bath transfer board. Four fixing points, no tools needed.
Measurements: 73x35 cm.



130 kg.



3'75 kg.



36x78x5 cm



2282

Bath transfer board with swivel foam seat. Four fixing points.
Quick assembly, no tools needed. It fits in baths from 43 to 68 cm. Size: 73x36cm. Seat diameter: 33cm.



130 kg.



3'75 kg.



38x76x8 cm



2198

Aluminium bath transfer bench, adjustable in height. Two legs are placed inside the bath, this way, allows the user to stay, get in and get out the bath without any efforts.



100 kg.



4'7 kg.



78x45x13 cm



AYUDAS
García 1880
MOVILIDAD

2296

Adjustable aluminium bath anti-slip seat. Seat size: 32x34.



110 kg.



37x35x9 cm



2297

U shaped aluminium bath anti-slip seat.



120 kg.



2294

Aluminium height adjustable shower bench with draining holes and anti-slip rubbers.



1'2 kg.



110 kg.



41x30x21 cm



2295

Aluminium adjustable stool to be used either in the shower or bathtub.

22931



Aluminium corner adjustable shower stool.



136 kg.



39x42x13 cm



No tools required.



110 kg.



1'45 kg.



42x37x8 cm



2203

Swivel bath chair with wheels.



53 cm.



56 cm.



55x55x40 cm



100 kg.



5'8 kg.



130 kg.



1'3 kg.



37x44x10 cm

2293

Aluminium corner shower stool.
Height adjustable: 43-53 cm. Seat
dimensions: 36x40 cm.

2299

Aluminium stool with soft swivel seat. Height
adjustable and suction rubbers.




100 kg.



2'1 kg.



39x42x13 cm



22951

Round stool with handles. Height
adjustable: 35.5-50.5 cm. Seat
diameter: 32cm. No tools needed.



136 kg.



1'3 kg.



34x48x7 cm



2201

Aluminium height adjustable swivel
chair with flip-up armrests and padded
backrest and seat. Seat with six
positions. No tools needed.



52 cm.



52 cm.



52x35x43 cm.



136 kg.



5'4 kg.



Bathroom aids for people with reduced mobility or disability can remain seated during the bath or shower.



2190

Aluminium height adjustable bath chair with backrest.



2191

Aluminium height adjustable bath stool.



2193

Aluminium height adjustable U shape bath chair with backrest.



2194

Aluminium height adjustable U shape bath stool.



21921 Epoxy painted steel foldable bath stool.



2192

Aluminium foldable bath stool. Height adjustable.



2196

Anodized aluminium height adjustable bath chair. Very light, with detachable armrests, backrest and draining holes.



46 cm.



53 cm.



64x12x40 cm



110 kg.



3,5 kg.



2197

Anodized aluminium height adjustable bath stool. Very light, with detachable armrests, backrest and draining holes.



46 cm.



53 cm.



36x64x13 cm



110 kg.



2,3 kg.



Asiento: 34x40 cm.

2290

Height adjustable bath chair with soft touch backrest and seat.



42 cm.



49 cm.




34x59x15 cm



100 kg.

2291

Height adjustable bath stool with soft touch seat.



2 kg.



2298

Aluminium foldable height adjustable bath chair.



56 cm.



52 cm.



52x78x21 cm



110 kg.



22981

Aluminium foldable bath chair with wheels, U shape ergonomic seat and flip up armrests and footrest. Wheels: 70 mm. Brakes on rear wheels.



53 cm.



60 cm.



107x51x25 cm.



90 kg.



6'0 kg.

21901

Deluxe aluminium shower chair with backrest and handles. Seat: 58x 36. Adjustable height. Max. Weight: 110 kg.



40 cm.



59 cm.



61x37x16 cm.



110 kg.



It is recommended that the installer values which screws to use when fixing the product. They will depend on the type of wall or the place where it is installed.

2188

Steel foldable wall mounted shower seat.



50x43x17 cm.



MAX.
110 kg.



5'1 kg.



21881

Deluxe wall mounted shower seat with support legs. Height adjustable: 46-53 cm.



40x46x10 cm.



MAX.
110 kg.



2185

Wall mounted U-shaped shower seat with support legs.



48x61x15 cm.



MAX.
136 kg.



2189

Deluxe wall mounted shower seat.



37x39x10 cm.



MAX.
110 kg.



Rubber mats for bathtub and shower

Made to prevent possible falls in areas where there is a risk of slipping, providing safety to the user.


2400

Non-slip bath mat. 35x75 cm.



2401

Non-slip bath mat. 40x104 cm.



Commode chairs



Daily aids often used by those who are unable to cover the distance to the bathroom due to mobility problems. They can also be placed directly on the toilet.

2170

Chromed steel height adjustable and foldable commode chair.



51 cm.



55 cm.



58x64x20 cm.



135 kg.



6'5 kg.



2171

Aluminium height adjustable foldable commode chair with wheels.



49 cm.



60 cm.



51x83x23 cm.



100 kg.



5'2 kg.



21082

Hi-Riser foldable aluminium walker with removable bucket.



60 cm.



50 cm.



64x11x75 cm.



110 kg.



2'5 kg.



2108

Hi-Riser foldable aluminium walker.

21081

Accessory bucket to be used as a toilet.



2211

Height adjustable commode chair. Steel frame. Upholstered and padded seat and backrest.



50 cm.



61 cm.



60x50x11 cm.
3 pieces



135 kg.



8'6 kg.



2212

Height adjustable commode chair. Steel frame. Upholstered and padded seat and backrest. Detachable frame. Easy assembly, no tools needed.



52 cm.



58 cm.



60x33x49 cm.



100 kg.



8'7 kg.



2209

Commode chair with wheels and detachable footrest. Frame made of steel. Padded seat. Rear wheels with brakes.



87 cm.



55 cm.



45x53x32 cm.



135 kg.



13 kg.





2230

Commode chair with wheels. Detachable footrest and flip up armrests. Frame made of steel. Padded seat. Rear wheels with brakes.



56 cm.



56 cm.



35x56x57 cm.



100 kg.



12 kg.



Shower commode chairs



Daily aids often used by those who are unable to cover the distance to the bathroom due to mobility problems. They can also be placed directly on the toilet.

2208

Height adjustable commode chair with wheels, flip up armrests and footrest.



73 cm.




56 cm.



95x60x23 cm.



150 kg.



2169

Aluminium commode chair with wheels and brakes on rear wheels.



22091

Aluminium commode chair with wheels and foldable footrest. Brakes on front and rear wheels.



21691

Aluminium commode chair with wheels and brakes on rear wheels.



22061

Aluminium tilting commode chair with flip-up armrest, height-adjustable footrest and adjustable padded headrest. Includes bucket and 4 wheels of 127 mm with brakes.



Quick Release Self-Propelled Shower Chairs



52 cm.



64 cm.



68x60x15 cm.
58x70x27 cm.



100 kg.



12'9 kg.



Código
SIRPO



Código
HERA-CV
SRM 040A-504



2234

Aluminium frame with PU backrest.

Quick release system for the rear wheels.

120 mm front wheels and 600 mm rear wheels. It also comes with 120 mm rear replaceable wheels.

Flip-up armrest. Adjustable and detachable footrests.

Seat width-depth: 45.5x41 cm. Total height: 104 cm. Seat-floor height: 53 cm.



107 cm.



69 cm.



70x61x40 cm.



150 kg.



18 kg.



22341

GRP structure (reinforced nylon with glass fibre).

Bucket included but can be used directly on the WC.

Quick release system for the rear wheels.

120mm front wheels and 600mm rear wheels.

It also comes with 120 mm rear replaceable wheels.

Folding armrest and removable footrest.

53 cm seat height. Total height 97cm.

Seat dimensions: 45x44 cm.

